

Robust Colony PCR from Multiple *E. coli* Strains using OneTaq[®] Quick-Load[®] Master Mixes

Yan Xu, Ph.D., New England Biolabs, Inc.

Introduction

Colony PCR is a commonly used method to quickly screen for plasmids containing a desired insert directly from bacterial colonies. This method eliminates the need to culture individual colonies and prepare plasmid DNA before analysis. However, the presence of bacterial cell contents and culture media in colony PCR reactions often results in polymerase inhibition. A robust polymerase is required to perform colony PCR with high efficiency in many different bacterial strains.

OneTaq DNA Polymerase, an optimized blend of *Taq* and Deep Vent_R[™] DNA polymerases, has been formulated for robust yields with minimal optimization. This robustness makes OneTaq ideal for use in demanding applications, such as colony PCR.

Furthermore, the OneTaq Quick-Load Master Mix product format increases the ease-of-use for colony PCR. The master mix formulation contains dNTPs, MgCl₂, buffer components and stabilizers, as well as two commonly used tracking dyes for DNA gels. On a 1% agarose gel in 1X TBE, Xylene Cyanol FF migrates at ~4 kb and Tartrazine migrates at ~10 bp. Both dyes are present in concentrations that do not mask any co-migrating DNA bands.

General Protocol

1. Transform ligation mix or other plasmid-containing reaction mixture into the desired bacterial strain, and incubate agar plates overnight at the appropriate temperature.
2. Set up 50 µl reactions as follows:

OneTaq Master Mix	25 µl
PCR primer	200 nM
H ₂ O	to 50 µl

Note: If OneTaq Hot Start Quick-Load 2X Master Mix is used, reactions can be set up at room temperature. If OneTaq Quick-Load 2X Master Mix is used, reactions should be set up on ice.

3. Use a sterile toothpick to pick up individual colonies and dip into each reaction tube.
4. As soon as the solution looks cloudy, remove the toothpick. To create a stock of each individual colony either:
 - a.) Dip the toothpick into 3 ml growth media with appropriate antibiotics and culture overnight.
 - or
 - b.) Streak the toothpick onto another agar plate containing the appropriate antibiotics and grow overnight.

(see other side)

DNA CLONING

DNA AMPLIFICATION & PCR

EPIGENETICS

RNA ANALYSIS

LIBRARY PREP FOR NEXT GEN SEQUENCING

PROTEIN EXPRESSION & ANALYSIS

CELLULAR ANALYSIS

Materials

- Well-isolated bacterial colonies, ideally 1-2 mm in diameter
- Sterile toothpicks or pipette tips
- Additional agar plate, or culture tubes with growth media for retention of original colonies.
- OneTaq Quick-Load 2X Master Mix with Standard Buffer (M0486) or OneTaq Hot Start Quick-Load 2X Master Mix with Standard Buffer (M0488)*
- Sterile H₂O
- PCR primers

**For amplicons with a GC content over 65% GC, OneTaq Quick-Load 2X Master Mix with GC Buffer or OneTaq Hot Start Quick-Load 2X Master Mix with GC Buffer may be used.*

5. Transfer reactions to a PCR cyclers, and perform PCR following the guidelines below for cycling conditions:

Initial denaturation:	
94°C	2 minutes
30 cycles:	
94°C	15–30 seconds
45–68°C	15–60 seconds
68°C	1 minute per kb
Final hold:	
68°C	5–10 minutes
10°C	hold


6. Load 4–6 µl of each PCR reaction directly onto an agarose gel, alongside an appropriate DNA ladder.

Results

Colony PCR was performed in 2 separate experiments using the protocol described above, with the following colonies:

- Colonies obtained from transformation of a plasmid with a 4.5 kb insert into 18 different *E. coli* strains. Amplification of the plasmid insert was achieved in each case. One *Taq* Quick-Load 2X Master Mix with Standard Buffer and One *Taq* Hot Start Quick-Load 2X Master Mix with Standard Buffer were used.

Colony PCR of a 4.5 kb insert using One *Taq* and One *Taq* Hot Start Quick-Load 2X Master Mixes with Standard Buffer and 18 different *E. coli* strains



Reactions were set up according to the protocol and analyzed by agarose electrophoresis. Marker M is the 1 kb DNA Ladder (NEB #N3232)

Lane	Name	NEB #	Lane	Name	NEB #	Lane	Name	NEB #
1	NEB 10-beta	C3019	7	Lemo21(DE3)	C2528	13	T7 Express lysY	C3010
2	NEB 5-alpha	C2987	8	NiCo21(DE3)	C2529	14	T7 Express	C2566
3	NEB 5-alpha F1q	C2992	9	NEB Express F	C3037	15	T7 Express Crystal	C3022
4	dam-/dcm-	C2925	10	NEB Express	C2523	16	SHuffle® Express	C3028
5	NEB Turbo	C2984	11	T7 Express F	C3016	17	SHuffle T7 Express lysY	C3030
6	BL21(DE3)	C2527	12	T7 Express lysY/F	C3013	18	SHuffle T7 Express	C3029

DNA CLONING

DNA AMPLIFICATION & PCR

EPIGENETICS

RNA ANALYSIS


LIBRARY PREP FOR NEXT GEN SEQUENCING

PROTEIN EXPRESSION & ANALYSIS

CELLULAR ANALYSIS

2. Colonies from *E. coli* library clones with inserts ranging from 0.8 kb to 10 kb. OneTaq Quick-Load 2X Master Mix with Standard Buffer was used, and results illustrate the robustness of the OneTaq Quick-Load 2X Master Mix in this application.

Colony PCR of library clones with inserts of 0.8 kb – 10 kb, using OneTaq Quick-Load 2X Master Mix



Ordering Information:

PRODUCT	NEB #	SIZE	PRICE
OneTaq® DNA Polymerase	M0480S/L/X	200/1,000/5,000 units	38,- € / 152,- € / 600,- €
OneTaq® 2X Master Mix with Standard Buffer	M0482S/L	100/500 rxns	38,- € / 152,- €
OneTaq® 2X Master Mix with GC Buffer	M0483S/L	100/500 rxns	38,- € / 152,- €
OneTaq® Quick-Load® 2X Master Mix with Standard Buffer	M0486S/L	100/500 rxns	39,- € / 156,- €
OneTaq® Quick-Load® 2X Master Mix with GC Buffer	M0487S/L	100/500 rxns	39,- € / 160,- €
OneTaq® Hot Start DNA Polymerase	M0481S/L/X	200/1,000/5,000 units	78,- € / 312,- € / 1248,- €
OneTaq® Hot Start 2X Master Mix with Standard Buffer	M0484S/L	100/500 rxns	65,- € / 268,- €
OneTaq® Hot Start 2X Master Mix with GC Buffer	M0485S/L	100/500 rxns	65,- € / 268,- €
OneTaq® Hot Start Quick-Load® 2X Master Mix with Standard Buffer	M0488S/L	100/500 rxns	68,- € / 278,- €
OneTaq® Hot Start Quick-Load® 2X Master Mix with GC Buffer	M0489S/L	100/500 rxns	68,- € / 278,- €
OneTaq® RT-PCR Kit	E5310S	30 rxns	149,- €

Featured in this Application Note.

Featured in this Application Note.

Summary

OneTaq and OneTaq Hot Start Quick-Load Master Mixes provide reliable performance in colony PCR, and are compatible with multiple *E. coli* strains. Reliable performance has been seen with amplicons up to 10 kb. The Quick-Load format offers additional convenience by enabling direct loading of the PCR reaction onto an agarose gel for analysis. Lastly, the Hot Start formulation provides additional functionality by reducing interference from primer-dimers and secondary amplification products.

NEW ENGLAND BIOLABS®, ONETAQ®, QUICK-LOAD® and SHUFFLE® are registered trademarks of New England Biolabs, Inc.

DEEP VENT™ is a trademark of New England Biolabs, Inc.



 NEW ENGLAND
BioLabs GmbH
enabling technologies in the life sciences

www.neb-online.de

New England Biolabs GmbH, Brüningstr. 50, Geb. B852, 65926 Frankfurt/Main, Germany

Tel: +49/(0)69/305-23140 Toll Free: (Germany) 0800/246-5227 Toll Free: (Austria) 00800/246-52277 Fax: +49/(0)69/305-23149 e-mail: info.de@neb.com