

ONE TAQ DNA POLYMERASE

ENZYMES & KITS FOR PCR




be INSPIRED
drive DISCOVERY
stay GENUINE

OneTaq DNA Polymerase –

The One You've been waiting for!

Choose OneTaq DNA Polymerase for all your amplification of standard, AT- and GC-rich templates, at a price-point that beats the competition. OneTaq DNA Polymerase is an optimized blend of Taq and Deep Vent DNA polymerases for use with routine and difficult PCR experiments. The 3'– 5' exonuclease activity of Deep Vent DNA Polymerase increases the fidelity and robust amplification of Taq DNA Polymerase. The OneTaq Reaction Buffers and High GC Enhancer have been formulated for robust yields with minimal optimization, regardless of a template's GC content (see additional data on page 4).

Achieve robust amplification for standard, AT- and GC-rich templates with OneTaq DNA Polymerase



Amplification of a selection of sequences with varying AT and GC content from human and *C. elegans* genomic DNA using OneTaq DNA Polymerase. GC content is indicated above gel. Marker M is the 1 kb DNA Ladder (NEB #N3232).

It's that easy to choose the right buffer for maximum performance:

OneTaq DNA Polymerase is supplied with two 5× buffers (Standard and GC), as well as a High GC Enhancer solution to ensure maximum performance for routine, AT- or GC-rich amplicons.

| Amplicon % GC content | Recommended default buffer | Optimization Notes |
|-----------------------|---------------------------------|--|
| <50% GC | OneTaq Standard Reaction Buffer | Adjust annealing temperature, primer/template concentration, etc., if needed. |
| 50–65% GC | OneTaq Standard Reaction Buffer | OneTaq GC Reaction Buffer can be used to enhance performance of difficult amplicons. |
| >65% GC | OneTaq GC Reaction Buffer | OneTaq GC Reaction Buffer with 10–20% OneTaq High GC Enhancer can be used to enhance performance of difficult amplicons. |



Advantages:

- Exceptional performance in endpoint PCR across a wide range of templates (AT-rich, Standard, GC-rich)
- Robust yields with minimal optimization
- Convenient product formats (stand-alone enzyme, master mixes, and Quick-Load formats)



Details & Applications:

Details

| | |
|----------------------------|-------------------------|
| Extension Rate | 1 kb/min |
| Amplicon Size | ≤ 6 kb |
| Fidelity | 2X Taq |
| Units/50 µl rxn | 1.25 units |
| Resulting Ends | 3' A/Blunt |
| 3'→5' Exonuclease Activity | Yes |
| 5'→3' Exonuclease Activity | Yes |
| Supplied Buffer | OneTaq Std Rxn Buffer |
| | OneTaq GC Rxn Buffer |
| Supplied Enhancer | OneTaq High GC Enhancer |

Product Formats

| | |
|------------------------------|-----|
| Hot Start Available | Yes |
| - Activation Required | No |
| Master Mix Available | Yes |
| Direct Gel-loading Available | Yes |
| PCR Kit Available | No |

Applications

| | |
|-------------------|-----|
| Routine PCR | Yes |
| SNP Detection | Yes |
| T/A, U/A Cloning | Yes |
| Colony PCR | Yes |
| High-Fidelity PCR | No |




For more information please visit www.neb.com/onetaq

OneTaq Quick-Load DNA Polymerase

– optimal for “Standard” PCR and direct gel-loading

For direct and fast agarose gel-loading after “standard” PCRs such as genotyping or colony-PCR etc., OneTaq DNA Polymerase is also available in a Quick-Load format. It is supplied with a density and tracking dye containing 5× OneTaq Quick-Load Reaction Buffer for direct gel loading in addition to the regular “colorless” 5× OneTaq Reaction Buffer.

Convenient direct gel-loading feature of OneTaq Quick-Load DNA Polymerase is not compromising on performance



Use your OneTaq PCR products in Sanger sequencing?

The dye in OneTaq Quick-Load 2X Master Mix Buffer doesn't interfere with Sanger sequencing. Prepare your samples with the fast and easy Exo-Cip Rapid PCR Cleanup Kit (#E1050) and proceed directly to Sanger sequencing.

Amplification of a variety of DNA targets demonstrates strong performance of the OneTaq Quick-Load DNA Polymerase. Product yield (bars, left axis) and purity (circles, right axis) were calculated via microfluidic analysis from triplicate reactions after 30 cycles of PCR. Standard deviation is indicated by error bars (yield) or shaded bands (purity). GoTaq was cycled according to manufacturer's recommendations.

Convenient Master Mix Formulations

– faster set-up, less pipetting, less errors & contaminations

OneTaq and OneTaq Hot Start DNA Polymerases are also available in convenient Master Mix and Quick-Load Master Mix formats. Master mix formulations include dNTPs, MgCl₂ and other buffers and stabilizers. The Quick-Load Master Mix formats allow for direct gel-loading. OneTaq Master Mix formulations are optimally suited for faster reaction set-up with less pipetting steps – they increase the reliability of each PCR and reduce the risk of contamination.

| OneTaq 2X Master Mixes | Art.-Nr. | convenient | gel-loading |
|--|------------|------------|-------------|
| with Standard Buffer | M0482S/L | • | |
| OneTaq Quick-Load 2× Master Mix with Standard Buffer | M0486S/L/X | • | • |

| OneTaq Hot Start 2X Master Mixes | Art.-Nr. | convenient | gel-loading |
|--|----------|------------|-------------|
| with Standard Buffer | M0484S/L | • | |
| with GC Buffer | M0485S/L | • | |
| OneTaq Hot Start Quick-Load 2× Master Mix with Standard Buffer | M0488S/L | • | • |
| OneTaq Hot Start Quick-Load 2× Master Mix with GC Buffer | M0489S/L | • | • |

Robust PCR results you can trust!

– unparalleled robustness on all tested templates!

While *OneTaq* DNA Polymerase works perfectly well on AT-rich DNA (see page 2), it's the challenging templates that unveil the true performance and quality of any PCR polymerase. Therefore, *OneTaq* DNA Polymerase and its Hot Start counterpart have been stringently and systematically tested under selected demanding conditions. In the figure below, we demonstrate some examples of the convincing performance of *OneTaq* on GC-rich templates.

Choose *OneTaq* DNA Polymerase for all your templates regardless of GC-content to benefit from the unrivaled robustness and reliability, so you can always trust your results!



Recommendation:

- Use *OneTaq* DNA Polymerase for all your templates
- Benefit from unrivaled reliability and performance
- Obtain PCR results you can really trust
- Choose the *OneTaq* format (e.g. convenient Master Mix) that suits you best
- Benefit from our low and fair prices

Better than the competition: A large, dark green dot represents the highest yield and purity.

| Polymerase | Additives | GC% | | | | | |
|---|----------------------|-----|----|----|----|----|----|
| | | 55 | 65 | 66 | 67 | 78 | 79 |
| OneTaq® DNA Polymerase (NEB) | None* | ● | ● | ● | ● | ● | ● |
| OneTaq HotStart DNA Polymerase (NEB) | None* | ● | ● | ● | ● | ● | ● |
| AmpliTaq Gold™ 360 DNA Polymerase (Thermo Fisher) | None | ● | · | ● | · | · | · |
| | 360 GC Enhancer | ● | · | ● | · | ● | ● |
| DreamTaq™ Hot Start DNA Polymerase (Thermo Fisher) | (Not provided) | ● | · | · | · | · | · |
| FastStart™ Taq DNA Polymerase (Roche) | None | ● | · | ● | · | · | · |
| | GC-RICH solution | ● | ● | ● | ● | ● | ● |
| GoTaq® G2 Hot Start Polymerase (Promega) | (Not provided) | ● | · | ● | · | ● | · |
| GoTaq Hot Start Polymerase (Promega) | (Not provided) | ● | · | ● | · | ● | · |
| HotStarTaq® DNA Polymerase (Qiagen) | Q-Solution | · | · | · | · | ● | ● |
| | None | ● | · | · | · | · | · |
| HotStarTaq Plus DNA Polymerase (Qiagen) | Q-Solution | · | · | · | · | ● | ● |
| | None | ● | · | · | · | · | · |
| iTaq™ DNA polymerase (Bio-Rad) | (Not provided) | ● | · | · | · | · | · |
| JumpStart™ Taq DNA Polymerase (Sigma) | None | ● | · | · | · | · | · |
| Platinum™ II Taq Hot-Start DNA Polymerase (Thermo Fisher) | Platinum GC Enhancer | ● | · | ● | ● | ● | ● |
| | None | ● | · | ● | · | · | · |
| Platinum Taq DNA Polymerase High Fidelity (Thermo Fisher) | (Not provided) | ● | · | · | · | · | · |
| Platinum Taq DNA Polymerase (Thermo Fisher) | None | ● | · | · | · | · | · |
| | KB Extender | · | ● | · | ● | ● | ● |
| Ex Taq DNA Polymerase, hot-start version (TaKaRa) | (Not provided) | ● | · | · | · | · | · |
| Titanium® Taq DNA Polymerase (TaKaRa) | (Not provided) | ● | · | · | · | · | · |

Yield (ng/ul)

· 0.0 ● 1.0 ● 2.0 ● 3.0 ● 4.0 ● ≥5.0

% Purity

0 100

* *OneTaq* products are supplied with *OneTaq* Standard Reaction Buffer and *OneTaq* GC Reaction Buffer. The GC reaction buffer was used to amplify the targets shown in the table. For other products, amplification reactions were conducted both with and without GC enhancers (if provided).

Amplification of a selection of high GC human genomic DNA targets demonstrates *OneTaq* performance. All polymerases were cycled according to manufacturer's recommendations, including the use of additives to enhance the amplification of targets with high GC content. Yield (dot size) and purity (color) of reaction product were quantified from triplicate reactions on a Perkin Elmer LabChip. A large, dark green dot represents the highest yield and purity.

OneTaq Hot Start DNA Polymerase

– room temperature reaction setup with no activation step

NEB's OneTaq Hot Start utilizes aptamer technology. This unique modified oligonucleotide binds to the polymerase through non-covalent interactions, blocking polymerase activity at temperatures below 45°C. The polymerase is activated during normal cycling conditions, allowing reactions to be set up at room temperature. OneTaq Hot Start DNA Polymerase does not require a separate high temperature incubation step to activate the enzyme. This ultimately shortens reaction times and increases ease of use.




Advantages:

- Allows room temperature reaction setup
- Does not require a separate activation step
- Compatible with standard Taq protocols

Recommended time for enzyme activation of commercially available Hot Start Taq products

| Manufacturer | Enzyme | Activation Step* | Hot Start Form |
|-----------------------------|-------------------|------------------|----------------|
| Applied BioSystems (Thermo) | AmpliTaq Gold 360 | 10', 95°C | Modified |
| Invitrogen (Thermo) | Platinum Taq | 30"-2', 94°C | Ab |
| Promega | GoTaq Hot Start | 2', 94-95°C | Ab |
| Qiagen | HotStarTaq | 15', 95°C | Modified |
| Roche | FastStart Taq | 4', 95°C | Modified |
| Sigma | JumpStart Taq | 1', 94°C | Ab |
| Thermo Fisher | Thermo-Start Taq | 15', 95°C | Modified |
| NEB | OneTaq | None | Aptamer |

* May include initial denaturation step



Extended room temperature incubation does not affect performance or specificity of OneTaq Hot Start DNA Polymerase

Amplification of a selection of sequences with varying GC content from human and *C.elegans* genomic DNA using OneTaq Hot Start DNA Polymerase. The presence or absence of an extended room temperature incubation does not affect performance. GC content is indicated above gel. Marker M is the 1 kb DNA Ladder (#N3232).



Get your free sample now!

Experience the advantages and request your **free sample** here!
www.neb.com/onetaq



Ordering information

| PRODUCTS | NEB # | SIZE |
|--|------------|---------------------------|
| OneTaq DNA Polymerase | M0480S/L/X | 200/1,000/ 5,000 units |
| OneTaq Quick-Load DNA Polymerase | M0509S/L/X | 100/500/ 2,500 units |
| OneTaq 2X Master Mix with Standard Buffer | M0482S/L | 100/500 rxns |
| OneTaq Quick-Load 2X Master Mix with Standard Buffer | M0486S/L/X | 100/500/ 2,500 rxns |
| OneTaq Hot Start DNA Polymerase | M0481S/L/X | 200/1,000/ 5,000 units |
| OneTaq Hot Start 2X Master Mix with Standard Buffer | M0484S/L | 100/500 rxns |
| OneTaq Hot Start 2X Master Mix with GC Buffer | M0485S/L | 100/500 rxns |
| OneTaq Hot Start Quick-Load 2X Master Mix with Standard Buffer | M0488S/L | 100/500 rxns |
| OneTaq Hot Start Quick-Load 2X Master Mix with GC Buffer | M0489S/L | 100/500 rxns |
| OneTaq RT-PCR Kit | E5310S | 30 rxns |
| OneTaq One-Step RT-PCR Kit | E5315S | 30 rxns |

Please ask for larger packing sizes or quantities: info.de@neb.com. Purchase of this product provides the purchaser with a non-exclusive license to use OneTaq DNA Polymerase products for research purposes only.

One or more of these products are covered by patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

AMPLITAQ GOLD, PLATINUM, THERMO-START, DREAMTAQ, are trademarks of Thermo Fisher Scientific.

FASTSTART is a trademark of Roche.

EX TAQ, TITANIUM are trademarks of TaKaRa Bio.

GOTAQ is a registered trademark of Promega.

HOTSTAR is a trademark of Qiagen.

ITAQ is a trademark of BIO-RAD.

JUMPSTART is a trademark of Sigma-Aldrich.



PCR Fidelity Estimator

For help with estimating the percentage of correct DNA copies, try PCRFidelityEstimator.neb.com



PCR Selector

For help with choosing the best polymerase for your PCR, try PCRselector.neb.com



Tm Calculator

For help with calculating annealing temperatures, choose our Tm Calculator at TmCalculator.neb.com

GERMANY & AUSTRIA

New England Biolabs GmbH
Brüningstr. 50, Geb B852
65926 Frankfurt/Main, Germany
Tel: +49/(0)69/305-23140
Fax: +49/(0)69/305-23149

Free Call: 0800/246 5227 (Germany)
Free Call: 00800/246 52277 (Austria)
info.de@neb.com

www.neb-online.de

FRANCE

New England Biolabs France SAS
Genopole Campus 1, Bâtiment 6
5 rue Henri Desbrüères
91030 Evry cedex, France
Tel.: 0800 100 632 (Customer Service)

Tel.: 0800 100 633 (Technical Service)
FAX.: 0800 100 610
info.fr@neb.com

www.neb-online.fr

HEADQUARTERS:

USA

New England Biolabs, Inc.
Telephone: (978) 927-5054

Toll Free (USA Orders): 1-800-632-5227

Toll Free (USA Tech): 1-800-632-7799

Fax: (978) 921-1350

info@neb.com

www.neb.com

www.neb.com



Distributors

For a complete list of all NEB subsidiaries and distributors, please go to:

www.neb.com/international-support

One or more of these products are covered by patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.